



# AGANI™ Hypodermic Needles

Smooth and easy puncture for patient comfort



# AGANI™ Hypodermic Needles

LATEX FREE

- Individually packaged in a sterile colour-coded blister pack for easy gauge identification.
- Suitable for all applications.
- Reassurance provided by visualisation of flashback through semi-transparent colour-coded needle hub. Thin wall needles allow you to choose a smaller gauge size with increased flow rates over regular walled needles.

## Agani™ Hypodermic Needles

Product Code	Colour	Wall Type	Gauge	Length	Carton
				mm	
AN*1838R1		Thin Wall	18G	38	10000
AN*1925R1		Thin Wall	19G	25	10000
AN*1938R1		Thin Wall	19G	38	10000
AN*2025R1		Thin Wall	20G	25	10000
AN*2038R1		Thin Wall	20G	38	10000
AN*2116R1		Thin Wall	21G	16	10000
AN*2125R1		Thin Wall	21G	25	10000
AN*2138R1		Thin Wall	21G	38	10000
AN*2219R1		Thin Wall	22G	19	10000
AN*2232R1		Thin Wall	22G	32	10000
AN*2238R1		Thin Wall	22G	38	10000
AN*2316R1		Thin Wall	23G	16	10000
AN*2325R1		Thin Wall	23G	25	10000
AN*2332R1		Thin Wall	23G	32	10000
AN*2425R1		Regular Wall	24G	25	10000
AN*2516R1		Thin Wall	25G	16	10000
AN*2525R1		Thin Wall	25G	25	10000
AN*2613R1		Regular Wall	26G	13	10000
AN*2623R1		Regular Wall	26G	23	10000
AN*2713R1		Regular Wall	27G	13	10000
AN*2716R1		Regular Wall	27G	16	10000
AN*2719R1		Regular Wall	27G	19	10000
AN*3013R1		Regular Wall	30G	13	10000

Shelf Pack = 100

## Needle Length: 50mm

Product Code	Colour	Wall Type	Gauge	Length	Carton
				mm	
AN*1850R1		Thin Wall	18G	50	5000
AN*1950R1		Thin Wall	19G	50	5000
AN*2150R1		Thin Wall	21G	50	5000
AN*2250R1		Thin Wall	22G	50	5000

Shelf Pack = 50

## Drawing Up Needle (Blunt Needle)

AN*1838N		Thin Wall	18G	38	10000
----------	--	-----------	-----	----	-------

Shelf Pack = 100

### Terumo Australia Pty Ltd

Macquarie Park NSW 2113 Australia  
Mount Albert Auckland 1025 New Zealand  
E: cs\_australia@terumo.co.jp

T: 1800 837 866  
T: 0800 66 77 57

AGANI™ trademark of Terumo Corporation.  
© 2020 Terumo Australia Pty Ltd. GHP200603D651V6 | PBD00349

